

X3^{plus}

X3Plus cart has a gas spring guided movement that makes height adjustments fast and easy.

Performance and comfort, without compromises.



Performance
and comfort,
beyond expectations

Top Accessories

Pole attachment
ø1,6cm



Bar code holder



PSU* holder



W 6, D 18, H 3,8 cm

Rail 1x2,5x30cm



Pole
H center 33,5 cm
Vesa 75/100
Tilting
Rotating



Double Pole***
H1 Center 45,5 cm
H2 Center 19 cm
Vesa 75/100
Tilting*
Rotating*
* upper only

Consoles

Console* in HPL laminate



W 54,5, D 40,5

Console* in coated steel



W 38, D 40cm

PC console
in Plastic with HPL top



W 51, D 40,5 H 8,5 cm

Note book console
in Plastic with HPL top (locker)




W 51, D 40,5 H 8,5 cm

Mouse plate Left/Right



W 15, D 17 cm

Qwerty drawer 1
W 45cm



Qwerty drawer 2
W 45cm



Column Accessories

Basket



W 35, D 25, H 10,5 cm


Humidifier bracket



Rail 1x2,5cm



wired Basket



W 38, D 25,5, H 14 cm

Drawer



W 35, D 25, H 10,5 cm



Silicone bellow



H adj. 25 cm

Oxygen support



3 drawers module
with dividers



W 45, D 36,5, H 43 cm

Tray



W 45, D 38 cm
W 35, D 28 cm


Under base accessories

PSU holder



W 6, D 18, H 3,8 cm

Brackets for Reomed 300/600W



Single wheel

ø 125 mm
standard
or ESD



Twin wheel

ø 125 mm
standard
or ESD



Technical specifications

Overall dimensions

540mm (W) x 580mm (D)
870mm (H)
without console

Packaging dimensions

680mm (W) x 774mm (D)
1100mm (H)

Weight

From 22 kg depending on the
configuration

CE Declaration of conformity

EN 60601-1
Directive 93/42/CE
medical device