


## X1






Handy, agile  
and complete,  
For all your needs

X1 Platform is a compact and safe equipment cart with a considerable loading capacity. The four lockable wheels of the base are responsible for an excellent manoeuvrability in narrow spaces and even in difficult conditions.








The aluminium extrusion column can be used for cabling and to connect different medical accessories, as drawer, shelf and basket that can be fixed at any position of the column.



## Top Accessories

- ☐ **Pole attachment**  
ø1,6cm  

- ☐ **Bar code holder**  

- ☐ **PSU\* holder**  
W 6, D 18, H 3,8 cm  

- ☐ **Rail 1x2,5x30cm**  

- ☐ **Pole**  
H center 33,5 cm  
Vesa 75/100  
Tilting  
Rotating  


## Consoles

- ☐ **Tablet support**  
10"  

- ☐ **Console**  
in HPL laminate, handle in ABS  
W 49, D 40,5 cm  

- ☐ **Console**  
in coated steel handle in ABS  
W 38, D 35,5 cm  

- ☐ **Qwerty drawer 1**  
W 45cm  

- ☐ **Console with VESA**  
75/100 Tilting  

- ☐ **Mouse plate Left/Right**  
W 15, D 17 cm  

- ☐ **Qwerty drawer 2**  
W 45cm  


## Column Accessories

- ☐ **Basket**  
  
W 35, D 25, H 10,5 cm
  - ☐ **Wired Basket**  
  
W 38, D 25,5, H 14 cm
  - ☐ **Drawer**  
  
W 35, D 25, H 10,5 cm

**Console Height 960cm**   
with ø100 wheels










**X1**  
**fixed**

- ☐ **Humidifier bracket**  

  - ☐ **Rail 1x2,5cm**  

  - ☐ **Tray**  
  
W 35, D 28 cm

## Under base accessories

- ☐ **PSU holder**  
W 6, D 18, H 3,8 cm  

- ☐ **Twin wheel**  
ø 75mm  

- ☐ **Twin wheel**  
ø 100mm  
standard  
or ESD  

- ☐ **Single wheel**  
ø 100mm  

- ☐ **Bumper**  
☐   
☐   
☐ 

## Technical specifications

**Overall dimensions**  
515mm (W) x 520mm (D)  
1010mm (H)

**Packaging dimensions**  
800mm (W) x 600mm (D)  
300mm (H)

**Weight**  
From 12 kg depending on the  
configuration

**CE Declaration of conformity**  
EN 60601-1  
Directive 93/42/CE