



**X1 Plus**  
Performance  
and ergonomy


**X1 Plus**, thanks to its gas spring adjustable height and its extended series of accessories, becomes the new standard for light medical trolleys. Ergonomy comes first: the new lifting column will enhance your performance and comfort in every situation.



## Top Accessories




## Consoles



## Column Accessories



## Under base accessories



## Technical specifications

**Overall dimensions**  
515 mm (W) x 520 mm (D)  
850/1070 mm (H)  
without console

**Packaging dimensions**  
800 mm (W) x 600 mm (D)  
300 mm (H)

**Weight**  
From 15 kg depending on the configuration

**CE Declaration of conformity**  
EN 60601-1  
Directive 93/42/CE  
medical device