

X3

A robust and
essential support
for your device

X3 platform stands out thanks to the care for design details and robustness. The console, with mechanical height adjustment from 70 to 110 cm and the aluminium column with a sliding accessories system, makes X3 a robust and essential support for your device.



Top Accessories * dimensions can be customized

☐ **Pole attachment**
ø1,6 cm



☐ **Bar code holder**



☐ **PSU* holder**



W 6, D 18, H 3,8 cm

☐ **Rail 1x2,5x30 cm**



☐ **Pole**

H center 33,5 cm
Vesa 75/100
Tilting
Rotating



☐

Double Pole

H1 Center 45,5 cm
H2 Center 19 cm
Vesa 75/100
Tilting*
Rotating*
* upper only



Consoles

☐ **Console*** in HPL laminate



W 54,5, D 40,5

☐ **Console*** in coated steel



W 38, D 40 cm

☐ **PC console**
in Plastic with HPL top



W 51, D 40,5 H 8,5 cm

☐ **Note book console**
in Plastic with HPL top (locker)



W 51, D 40,5 H 8,5 cm

☐ **Mouse plate** Left and Right



W 51, D 17 cm

☐ **Mouse plate** Left/Right




W 15, D 17 cm

☐ **Qwerty drawer 1**
W 45 cm



☐ **Qwerty drawer 2**
W 45 cm



Column Accessories

☐ **Basket**



W 35, D 25, H 10,5 cm

☐ **Wired Basket**



W 38, D 25,5, H 14 cm

☐ **Drawer**



W 35, D 25, H 10,5 cm

☐ **Humidifier bracket**



☐ **Rail 1x2,5 cm**



Console height
from 70 cm to 110 cm



X3

Lateral sliding

☐ **Oxygen support**



☐ **3 drawers module**
with dividers



W 45, D 36,5, H 43 cm

☐ **Tray**



W 45, D 38 cm
W 35, D 28 cm


Under base accessories

☐ **PSU holder**



W 6, D 18, H 3,8 cm

☐ **Brackets** for Reomed 300/600W



☐ **Single wheel**
ø125 mm
standard
or ESD



☐ **Twin wheel**
ø125 mm
standard
or ESD



Handle and bumpers



Technical specifications

Overall dimensions

540mm (W) x 580mm (D) 890 (H)
without console

Packaging dimensions

800mm (W) x 600mm (D)
300mm (H)

Weight

From 20 kg depending on the
configuration

CE Declaration of conformity

EN 60601-1
Directive 93/42/CE
medical device