



MULTIX

Tailored Medical Trolley

Endless Possibilities

TO CUSTOMIZE YOUR TROLLEY



MULTIX is designed to expand the possibilities of your devices. It provides a broad spectrum of solutions: from a range of ready-to-use accessories to a fully customised set of tools developed for your own trolley.

MULTIX provides smart solutions specifically designed to help your equipment to face new challenges, improve your performances and achieve better results.

X1



X1Plus



X3



X3Plus





Hygiene
Cleanable surfaces



Design
Safe, intuitive and essential



Tailored
We develop custom solutions for every need



Ergonomy
User-centered design:
Comfortable for different users



Stability
Optimized footprint design



Multifunctional
Different configurations
for different needs



Robustness
Aluminium and steel structure




Smart and Agile: **X1 Platform**

X1 Platform is a compact and safe equipment cart with a considerable loading capacity. The four lockable wheels of the base ensure an excellent manoeuvrability in narrow spaces and even in difficult conditions.

The extruded aluminium column can be used for cabling and to connect different medical accessories, such as drawer, shelf and basket that can be fastened in any position of the column.

X1

Handy, agile and complete,
for all your needs



MULTIX

Tailored Medical Trolley



Compact Footprint
515x515 mm



Example of configuration




NEW




Example of configuration
for **tablet**

X1 Plus

Performance
and ergonomy



MULTIX
Tailored Medical Trolley





Hygiene
Cleanable surfaces



Design
Safe, intuitive and essential



Tailored
We develop custom solutions
for every need



Ergonomy
User-centered design:
comfortable for different users




Stability
Optimized footprint design



Multifunctional
Different configurations
for different needs



Robustness
Aluminium and steel structure




Turn the trolley **upside-down**
and appreciate the care for
the hidden details.

Strong and ergonomic: **X3 Platform**


X3 platform stands out thanks to the care for design details.
X3 platform will have a significant impact on your workflow
thanks to its robustness.

X3

A robust and essential support for your device.



Sliding height adjustment, perfect for installing your device in the most ergonomic position.



Aluminium column 70x140 mm with a sliding accessories system

MULTIX
Tailored Medical Trolley



Example of configuration for device




Example of configuration with **HPL laminated console**

ø125 mm
Base with single wheels.
Twin wheels optional.

X3 Plus

Performance and comfort,
beyond expectations



Customized logo



Silicone bumpers

MULTIX

Tailored Medical Trolley



ø125 mm
Choice between single or twin castors, conductive or standard





MULTIX provides dedicated solutions for any need: the notebook console gives an ergonomic workspace for notebook up to 17" with an easy lock system.



Your tailored trolley

READYMADE
ACCESSORIES

MULTIX
Tailored Medical Trolley

X1, X1^{Plus} Consoles

Tablet support
10"



Console *
in HPL laminate, handle in ABS




W 49, D 40,5 cm *

Console with VESA
75/100 Tilting



Console *
in coated steel handle in ABS



W 38, D 35,5 cm

Top Accessories

PSU* holder



W 6, D 18, H 3,8 cm

Rail 1x2,5x30 cm




Pole attachment
ø1,6 cm



Bar code holder



Qwerty drawer 1
W 45 cm



Qwerty drawer 2
W 45 cm



PC console
in Plastic with HPL top



W 51, D 40,5 H 8,5 cm

Note book console
in Plastic with HPL top (locker)



W 51, D 40,5 H 8,5 cm

Mouse plate Left/Right



Mouse plate Left/Right



W xx, D xx cm



Double Pole***

H1 Center 45,5 cm

H2 Center 19 cm

Vesa 75/100

Tilting*

Rotating*

* upper only



Pole

H center 33,5 cm

Vesa 75/100

Tilting

Rotating

Choose your color accent



X3, X3^{Plus} Consoles

Console* in HPL laminate



W 54,5, D 40,5


Console* in coated steel



W 38, D 40 cm

X1, X1^{Plus} Platform

X1 platform



X1
fixed



X1^{Plus}
Gas spring assistant
H adj. 22 cm


Column Accessories

Basket



W 35, D 25, H 10,5 cm

wired Basket



W 38, D 25,5, H 14 cm

Drawer



W 35, D 25, H 10,5 cm

Humidifier bracket



Rail 1x2,5 cm



3 drawers module***
with dividers



W 45, D 36,5, H 43 cm

Tray



W 45, D 38 cm*

W 35, D 28 cm

Oxygen support***




X3
lateral sliding

Silicone Bellow***
H. adj. 25 cm



X3^{Plus}
Gas spring assistant
H. adj. 29 cm

Under base accessories



Twin wheel**
ø 75 mm



Twin wheel
ø 100 mm
standard
or ESD



Single wheel
ø 100 mm

PSU* holder



W 6, D 18, H 3,8 cm

Brackets for Reomed 300/600W***



Single wheel***
ø 125 mm
standard
or ESD



Twin wheel***
ø 125 mm
standard
or ESD



MULTIX

Tailored Medical Trolley

We provide a complete design service that supports you
and your customers, not just a product.

MULTIX helps your equipment to improve your performances
and achieve better results because we offer a device-centered
and client-centered design service.

design by:

rezzonico|*design*



reddot design award



INTERNATIONAL DESIGN
EXCELLENCE AWARDS



Premio Compasso d'Oro ADI



A DESIGN AWARD
PLATINUM

For more than 25 years the Rezzonico Team has developed competitive
cutting-edge products and innovative design solutions
in the healthcare industry.

Multix Solutions s.r.l.
phone +39 031 4100263
fax +39 031 269994
info@multix-trolley.com
www.multix-trolley.com

