



ACCESSORIES



MULTI 

all the accessories you need

X1, X1T, X2, Platforms

endless


configurations

Console Accessories

2 or 3 accessories on each side can be added to the console.


The optional Qwerty and left/right mouse slide in underneath the front of the console.

**Customized accessories
can be designed on request!**






IV POLE




PROBE HOLDERS



BAR CODE HOLDERS



**QWERTY support
with LEFT/RIGHT Mouse**



Multix, Your Partner
in **Healthcare Excellence**


QWERTY SUPPORT

(W) 490mm
(D) 155mm

Max load 1Kg

ART. CODE:
025MLTACC

Available for:
X1, X1T Platforms




PIVOT HINGE FOR HEIGHT ADJ.VESA

Support to be mounted
on VESA +-90°

ART. CODE:
019MLTACC

Available for:
X1, X1T, Platforms





L-R MOUSEPAD


REQUIRES QWERTY SUPPORT (025MLTACC)

(W) 140 mm
(D) 170 mm

Max load 1Kg

ART. CODE:
017MLTACC

Available for:
X1, X1T Platforms




SIDE RAIL

10x25mm DIN rail

Length 30 cm
Max load 4Kg

ART. CODE:
028MLTACC

Available for:
Console H,L,M
All Platforms



CAMERA HOLDER

**Support to be mounted
on VESA with
Height adjustment
from 19" to 27"**

Camera Plate 150x150 mm

Max load 1,5Kg

ART. CODE:
010MLTACC

Available for:
All Platforms



SIDE RAIL PLASTIC CONSOLE

10x25mm DIN rail

Length 30 cm
Max load 3Kg

ART. CODE:
051MLTACC

Available for:
X1, X1T Platforms
Console C, CE, N, NL, CL



IV POLE

(H) 350mm
Max load 4Kg
2Kg per hook

ART. CODE:
015MLTACC

Available for:
Console H,L,M
All Platforms




SANITIZER GEL HOLDER

Inner dimensions:
80x80mm
Height 112mm
Max load 1Kg

ART. CODE:
013MLTACC

Available for:
All Platforms




PROBE HOLDER C-D-E

PMMA Transparent

Max load 3Kg

Available for:
All Platforms



C

ART. CODE:
033MLTACC



D

ART. CODE:
034MLTACC



E

ART. CODE:
035MLTACC


PROBE HOLDER G

Stell Rubbered

Max load 3Kg

ART. CODE:
049MLTACC

Available for:
All platforms




PROBE HOLDER A

Inner ø min 16 mm
Inner ø max 55 mm
Inner Height 42mm

Max load 1Kg

ART. CODE:
023MLTACC


Available for:
All platforms



POLE ATTACHMENT A

Max load 4Kg

Available for:
Console H,L,M
All Platforms



ø 18mm

ART. CODE:
039MLTACC

ø 16mm

ART. CODE:
040MLTACC

ø 13mm

ART. CODE:
041MLTACC

POLE ATTACHMENT B

Max load 4Kg

Available for:
X1, X1T Platforms
Console C, CE, N,
NL, CL



ø 18mm

ART. CODE:
058MLTACC

ø 16mm

ART. CODE:
059MLTACC

ø 13mm

ART. CODE:
060MLTACC

PROBE HOLDER B

Inner ø min 20 mm
Inner ø max 55 mm
Inner Height 82mm

Max load 3Kg

ART. CODE:
024MLTACC

Available for:
All platforms






TABLET SUPPORT

Support to be mounted
on VESA

Max load 4Kg

ART. CODE:
066MLTACC

Available for:
All VESA 75/100
attachments




BARCODE TPM

Max load 2Kg

ART. CODE:
070MLTACC

Available for:
All Platforms




SCANNER MOUNTING PLATE

CR2500, CR2600, CR2700, ZEBRA CR8178,
DATALOGIC 4400, HONEYWELL 1902,1952,1962

Max load 2Kg

ART. CODE:
069MLTACC

Available for:
All Platforms




BARCODE A

Max load 2Kg

ART. CODE:
005MLTACC

Available for:
All Platforms




BARCODE B


Max load 2Kg

ART. CODE:
006MLTACC

Available for:
All Platforms




Aluminium column
with front/back
Sliding accessories system



All accessories
can slide along
the column and
securely lock
into place using a
pressure-based
system

Column
Accessories



Column Accessories




	DOUBLE	SINGLE
X1	2 MAX	3 MAX
X1 T	1 MAX	2 MAX
X1 NB	1 MAX	2 MAX

How many drawers can be mounted on each model?

1 RAIL
DIN 10x25mm



2 PC BASKET



3 PC DRAWER



SINGLE DRAWER



DOUBLE DRAWER



SINGLE DRAWER*

with KEY LOCKING SYSTEM

ART. CODE:
003MLTACC


Available for:



Inner dim.
(W)347mm
(D) 247mm
(H) 110 mm

Outer dim.
(W)465 mm
(D) 360 mm
(H) 134mm

Max load 3Kg



open drawer
with Optional dividers



compatible with SALTO or RFDI locker

DOUBLE DRAWER*

with KEY LOCKING SYSTEM

ART. CODE:
004MLTACC

Available for:
X1, X1T
Platforms




Outer dim.
(W) 465 mm
(D) 360 mm
(H) 268mm

Max load 6Kg



compatible with SALTO or RFDI locker

removable tray
ISO 40x30x10




*SINGLE AND DOUBLE DRAWERS
ARE AVAILABLE WITH OPTIONS FOR DIVIDERS OR ISO BASKETS

SINGLE PLASTIC DRAWER*

with KEY LOCKING SYSTEM

ART. CODE:
075MLTACC


Available for:
X1, X1T
Platforms



Inner dim.
(W) 360mm
(D) 280mm
(H) 106 mm

Outer dim.
(W) 434 mm
(D) 390 mm
(H) 156mm

Max load 3Kg



compatible with SALTO or RFDI locker

**open drawer
without Optional dividers**

PC DRAWER

Polycarbonate, finish: frosted

ART. CODE:
014MLTACC

Available for:
X1, X1T
Platforms




(W) 333 mm
(D) 243mm
(H) 101mm
Max load 3Kg

PC BASKET

Polycarbonate, finish: frosted

ART. CODE:
007MLTACC

Available for:
X1, X1T
Platforms



(W) 333mm
(D) 243mm
(H) 101mm
Max load 3Kg


WIRED BIN

(W) 315mm
(D) 100mm
(H) 135mm

Max load 3Kg

ART. CODE:
031MLTACC

Available for:
All Platforms




WIRED BASKET

(W) 380mm
(D) 280mm
(H) 142 mm

Max load 3Kg

ART. CODE:
053MLTACC

Available for:
X1, X1T Platforms



SHELF

Internal
dimensions
(W) 390mm
(D) 310mm

Max load 15Kg

Available for:
X1, X1T Platforms

ART. CODE:
032MLTACC



WIRED BASKET WITH CABLE HOLDER

(W) 370mm
(D) 252mm
(H) 119 mm

Max load 3Kg

ART. CODE:
054MLTACC

Available for:
All Platforms



SOCKET/PSU SUPPORT

(W) 96 mm
(H) 238 mm

Max load 3Kg

ART. CODE:
029MLTACC

Available for:
All Platforms



MEDICAL MULTI SOCKET

N.4 C14 DIN multi sockets

Main Power C13 DIN and Switch

(W) 121mm
(H) 370mm
(D) 100mm

ART. CODE:
073MLTACC

Available for:
All Platforms



PRINTER TRAY A

Inner dim.
(W) 450 mm,
(D) 380mm,

Max load 10Kg

ART. CODE:
020MLTACC

Available for:
X1 Platform



PC SUPPORT

Internal dimensions
(W) 440mm
(D) 185mm

Max load 10Kg

ART. CODE:
030MLTACC

Available for:
All Platforms



PEDAL SWITCH SUPPORT

Internal dim.
163 mm (W) × 41 mm (D) x 156 mm (H)

Max load 3Kg

ART. CODE:
048MLTACC

Available for:
All Platforms



COLUMN RAIL

(W) 300mm,

Max load 10Kg

ART. CODE:
055MLTACC

Available for:
All Platforms



DC POWER BATTERY SUPPORT

ready for QPD DC Box
DCPower4C battery

Housing
(W) 260 mm
(D) 105 mm
(H) 351mm

Max load 9Kg

ART. CODE:
012MLTACC

Available for:
All Platforms



HUMIDIFIER SUPPORT

Max load 3Kg

ART. CODE:
079MLTACC

Available for:
All Platforms



CABLE CHANNEL 47/70/77/84/110

Cable Channel 3 ways-470/720mm
About 6/8 cables

ART. CODE:
008MLTACC_47
008MLTACC_70
008MLTACC_77
008MLTACC_84
008MLTACC_110

Available for:
All Platforms




VESA COLUMN ATTACHMENT

VESA 75-100

Max load 6Kg

ART. CODE:
042MLTACC

Available for:
All Platforms




SINGLE CYLINDER HOLDER 3L-5L-10L

Max Diam. ø140mm

Max load 7Kg

ART. CODE:
038MLTACC

Available for:
X1 and X2 Platforms



Base/column Accessories

The following accessories are to be mounted either solely on the base or partially on both the base and the column.

Multix, Your Partner
in **Healthcare Excellence**

DOUBLE 5L CYLINDER SUPPORT

Max Diam. ø140mm

Max load 14Kg

ART. CODE:
072MLTACC

Available for:
X1 Platform




5L CYLINDER SUPPORT

Max Diam. ø140mm

Max load 7Kg

ART. CODE:
052MLTACC

Available for:
X1 T




DOUBLE 10L CYLINDER SUPPORT

Max Diam. $\varnothing 190\text{mm}$

Max load 20+20Kg

ART. CODE:
076MLTACC

Available for:
X2 Platform




DOUBLE 5L CYLINDER SUPPORT

Max Diam. $\varnothing 100\text{mm}$

Max load 5+5Kg

ART. CODE:
077MLTACC

Available for:
X2 Platform



PRINTER TRAY X1T

Inner dim.
(W) 450 mm,
(D) 380mm,

Max load 10Kg

ART. CODE:
021MLTACC

Available for:
X1T



PSU SUPPORT


UNDER BASE

Inner dim.
92 mm (W) $\times$ 232 mm (D) $\times$ 55 mm (D)

Max load 5Kg

ART. CODE:
064MLTACC

Available for:
X1, X1T Platforms



INSULATOR TRANSF. SUPPORT


IMEDe 300 3rd, IMEDe 600 3rd, IMEDe 600

Inner dim.
(H) 104 mm

Max load 9Kg

ART. CODE:
078MLTACC

Available for:
X1, X1T Platforms



INSULATOR TRANSF. SUPPORT


REOMED I 300 - REOMED I 600

Inner dim.
(H) 97,5mm

Max load 12Kg

ART. CODE:
036MLTACC

Available for:
X1, X1T Platforms





MLTCAT_ACC_01 (EN)



Scan the brochure



TOUCH POINT GROUP

Multix Solutions s.r.l.
phone +39 031 4100263

info@multix-trolley.com
www.multix-trolley.com

Via Padre Rovagnati, 2
22066 Mariano Comense
COMO - ITALY

Made in Italy



ACCESSORIES
CATALOGUE