



X1 Platform

PURE ALUMINIUM

MULTIX

X1 PLATFORM

Compact and Robust

Efficient and **Compact Design**

Our platform's compact design ensures it fits seamlessly in various medical settings, from busy hospital wards to compact clinics. Despite its small footprint, it boosts efficiency and keeps essential tools within reach, proving indispensable for healthcare professionals.

Durable and **Robust Construction**


Constructed from high-quality materials, the Multix platform is built to endure the rigors of healthcare environments, offering unwavering reliability and a long-lasting solution for any medical facility.

Adaptable to **Your Needs**

Recognizing the diversity of medical practices, **our platform offers versatile configurations**, detailed in the following pages. Beyond these, **we also provide customization** options to meet the specific demands of your facility, from adjustable compartments to unique add-on features.

The Multix platform is more than equipment—it's a partner in enhancing healthcare delivery. Discover how our compact, robust, and customizable solution can revolutionize your practice.

Multix: Your Partner in Healthcare Excellence



X1BS77M accessorized with:
SINGLE PLASTIC drawer (075MLTACC)
Column Rail (011MLTACC)

Swivelling $\pm 30^\circ$

Tilting $0^\circ + 40^\circ$

**Monitor or
Panel PC**

Up to 24"

Max Load: 12kg



X1BS84PF accessorized with:
PC Drawer (014MLTACC)
Shelf (032MLTACC)

Let us inspire you

Swivelling $\pm 30^\circ$

Tilting $0^\circ + 40^\circ$

**Monitor or
Panel PC**

Up to 24"
Max Load: 12kg



X1BS77HF accessorized with:

Basket	(007MLTACC)
Gel Holder	(013MLTACC)
Barcode B	(006MLTACC)

UP to 13 cm

Tilting $-15^\circ + 25^\circ$

**Monitor or
Panel PC**

Up to 24"
Max Load: 12kg



X1GT84CR accessorized with:

PC Drawer	(014MLTACC)
Qwerty Tray	(025MLTACC)
Mouse Pad	(017MLTACC)




Console for Notebook
UP to 17"



X1BT77N accessorized with:
ISO Drawer (004MLTACC)


Support VESA for tablet
UP to 13"



X1BS84T accessorized with:
Wired bin (031MLTACC)
SINGLE PLASTIC drawer (075MLTACC)

Let us inspire you


X1 Platform




70
77
84

X1N
Without
Counterweights

Accent Color
● **G** Gray
● **B** Blue





F	Fixed support VESA 75/100 * ø60 mm Pole, ** Oval
R	Height adjustable support
N	Notebook Enclosure
C	PC case
CE	PC case
L	HPL work shelf
H	HPL work shelf with handles
M	Metal plate
P	Without console
T	Tablet support
84	COLUMN (H 838mm) FOOT PRINT (W)520x(D)520 mm
70	COLUMN (H 694 mm) FOOT PRINT (W)520x(D)520 mm
77	COLUMN (H 765 mm) FOOT PRINT (W)520x(D)520 mm
T	TWIN WHEELS Ø 100 mm 4 with Brakes, 2 ESD
S	SINGLE WHEELS Ø 100 mm 4 with Brakes, 2 ESD

	Color	Wheels	Column	Console	Monitor Supp.
X1	G B	S T	84 77 70	C CE L H M P T	/ F R
	☐	☐	☐	☐	☐
X1N	G B	S T	84 77 70	N	
	☐	☐	☐	☐	



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Case (S09) inner size	(W)440,(D)280,(H)74 mm
Top Material	HPL with plastic case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70C-F	Blue	Single 100mm	700mm	995 mm	1351 mm	800x570x370	22Kg
X1-BS-77C-F	Blue	Single 100mm	770mm	1065 mm	1421 mm	800x570x370	22Kg
X1-BS-84C-F	Blue	Single 100mm	840mm	1135 mm	1491 mm	880x600x370	22Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Case (S09) inner size	(W)440,(D)280,(H)74 mm
Top Material	HPL with plastic case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70C-R	Blue	Single 100mm	770mm	995 mm	1396 mm	800x570x370	23Kg
X1-BS-77C-R	Blue	Single 100mm	770mm	1065 mm	1466 mm	800x570x370	23Kg
X1-BS-84C-R	Blue	Single 100mm	770mm	1065 mm	1536 mm	880x600x370	23Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY

X1 Platform



General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Case (S09) inner size	(W)440,(D)280,(H)74 mm
Top Material	HPL with plastic case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

MAX. RATING
10A, 250V A.C.
(15A, 250V A.C. UL&CSA)



Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70CE-F	Blue	Single 100mm	700mm	995mm	1351mm	800x570x370	22Kg
X1-BS-77CE-F	Blue	Single 100mm	770mm	1065mm	1421mm	800x570x370	22Kg
X1-BS-84CE-F	Blue	Single 100mm	840mm	1135mm	1491mm	880x600x370	22Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform



General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Case (S09) inner size	(W)440,(D)280,(H)74 mm
Top Material	HPL with plastic case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa height adj	130 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°

MAX. RATING
10A, 250V A.C.
(15A, 250V A.C. UL&CSA)



Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70CE-R	Blue	Single 100mm	770mm	995 mm	1396 mm	800x570x370	23Kg
X1-BS-77CE-R	Blue	Single 100mm	770mm	1065 mm	1466 mm	800x570x370	23Kg
X1-BS-84CE-R	Blue	Single 100mm	770mm	1065 mm	1536 mm	880x600x370	23Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

WWww

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70L-F	Blue	Single 100mm	700mm	905 mm	1230mm	800x570x370	20Kg
X1-BS-77L-F	Blue	Single 100mm	770mm	975 mm	1300mm	800x570x370	20Kg
X1-BS-84L-F	Blue	Single 100mm	840mm	1045 mm	1370 mm	880x600x370	20Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with front handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa height adj	130 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70L-R	Blue	Single 100mm	700mm	905mm	1306mm	800x570x370	21Kg
X1-BS-77L-R	Blue	Single 100mm	770mm	975mm	1376mm	800x570x370	21Kg
X1-BS-84L-R	Blue	Single 100mm	840mm	1045mm	1446mm	880x600x370	21 Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70H-F	Blue	Single 100mm	700mm	905 mm	1230mm	800x570x370	20Kg
X1-BS-77H-F	Blue	Single 100mm	770mm	975 mm	1300mm	800x570x370	20Kg
X1-BS-84H-F	Blue	Single 100mm	840mm	1045 mm	1370 mm	880x600x370	20Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

  *AVAILABLE ALSO IN GREY



X1 Platform


General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with front handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Vesa height adj	130 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°

Ref.	Accent Color	Wheel	Column Length	Console Height	C.VESA MAX Height	Packaging Dimensions mm	Weight
X1-BS-70H-R	Blue	Single 100mm	700mm	905mm	1306 mm	800x570x370	21Kg
X1-BS-77H-R	Blue	Single 100mm	770mm	975 mm	1376 mm	800x570x370	21Kg
X1-BS-84H-R	Blue	Single 100mm	840mm	1045 mm	1446 mm	880x600x370	21 Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

  *AVAILABLE ALSO IN GREY

X1 Platform



General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	440(W)x400(D) mm
Top Material	Steel
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg


Ref.	Accent Color	Wheel	Column Lenght	Console Height	Packaging Dimensions mm	Weight
X1-BS-70M	Blue	Single 100mm	700mm	900 mm	800x570x370	20Kg
X1-BS-77M	Blue	Single 100mm	770mm	970 mm	800x570x370	20Kg
X1-BS-84M	Blue	Single 100mm	840mm	1040 mm	880x570x370	20Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY

X1 Platform



General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Case (S09) inner size	(W)440,(D)280,(H)74 mm
Top Material	HPL with plastic case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg




Ref.	Accent Color	Wheel	Column Lenght	Console Height	Packaging Dimensions mm	Weight
X1N-BS-70N	Blue	Single 100mm	700mm	997 mm	800x600x320	15 Kg
X1N-BS-77N	Blue	Single 100mm	770mm	1067 mm	800x600x320	15 Kg
X1N-BS-84N	Blue	Single 100mm	840mm	1137 mm	800x600x320	15 Kg

(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Tablet vesa attachment	75/100
Compliant with	MDR 2017/745 EN 60601-1
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°

Ref.	Accent Color	Wheel	Column Lenght	Console Height	Packaging Dimensions mm	Weight
X1-BS-84P-F	Blue	Single 100mm	840mm	1375 mm	880x600x370	16 Kg


(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



*AVAILABLE ALSO IN GREY



X1 Platform



General information	
Base foot print	520x520mm
Max Load	50kg
Tablet vesa attachment	75/100
Compliant with	MDR 2017/745 EN 60601-1
Center Vesa	xxxx
Vesa max load	12kg
Vesa tilting	-0°+40°

Ref.	Accent Color	Wheel	Column Lenght	Console Height	Packaging Dimensions mm	Weight
X1-BS-84T	Blue	Single 100mm	840mm	1121 mm	880x600x370	16 Kg
X1N-BS-84T	Blue	Single 100mm	840mm	1121 mm	880x600x370	16 Kg


(T instead of S) When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

  *AVAILABLE ALSO IN GREY

Some examples of custom Configurations


Custom accessories for medical trolleys are essential for optimizing efficiency and organization in healthcare settings. These accessories, tailored to specific needs, include drawer dividers, lockable compartments, adjustable shelving, power outlets, waste bins, accessory hooks, documentation holders, and customized labeling. By integrating these features, medical facilities can streamline workflow, enhance patient care delivery, and ensure a clean and organized environment.

Multix, Your Partner
in **Healthcare Excellence**





Customized console
for dental device



**Customized
Switch Pedal
Support**




Custom solutions



Customized column
for NICU

**Customized console
and shelves**
for TECAR device.





Customized
for endoscopy.



Custom solutions

MULTIX

MLTCAT X1_03 (EN)



TOUCH POINT GROUP

Multix Solutions s.r.l.
phone +39 031 4100263

info@multix-trolley.com
www.multix-trolley.com

Via Padre Rovagnati, 2
22066 Mariano Comense
COMO - ITALY

Made in Italy



X1 Platform
CATALOGUE