




X1T Platform

NEXT LEVEL MEDICAL TROLLEY

MULTIX

X1T PLATFORM

Ergonomic and Agile

Our platform, the X1T series, epitomizes ergonomic excellence and agility in medical equipment. Thanks to its gas spring column that affords an impressive 35 cm of smooth vertical travel. This feature not only enhances user comfort but also facilitates swift manoeuvrability in dynamic medical environments.

Tailored to Your Application

Recognizing the diverse needs of medical practitioners, our X1T platform offers versatile configurations meticulously outlined in accompanying documentation. Moreover, we provide customizable options to meet your requirements, whether dimensional customization or specialized add-on features.

Experience the Difference

More than just equipment, the X1T platform is a reflection of our commitment to enhancing healthcare delivery. Discover how its ergonomic design, agility, and customizable features can revolutionize your practice, setting a new standard for efficiency and mobility in medical care.

Multix: Your Partner in Healthcare Excellence



X1T-BS-20H-R

Without accessories




X1TGT40CF accessorized with:

PC Drawer	(Code 014MLTACC)
Qwerty Tray	(Code 025MLTACC)
Mouse Pad	(Code 017MLTACC)
Probe Holder	(Code 023MLTACC)
Side Rail	(Code 028MLTACC)

Let us inspire you

areas of application


- Information Technology
- Telemedicine
- Endoscopy and Operating room
- ECG
- Ultrasound
- Rehabilitation




X1TBS20M accessorized with:

Wired Bin (Code 031MLTACC)

Printer Tray (Code 021MLTACC)



X1TBT30N accessorized with:
N. 02 AL Drawer (Code 003MLTACC)



X1Tall accessorized with:
N.01 ISO Drawer (Code 004MLTACC)
Qwerty Tray (Code 025MLTACC)
Mouse Pad (Code 017MLTACC)
Barcode B (Code 006MLTACC)
Side Rail (Code 028MLTACC)

Let us inspire you


X1T Platform




Accent Color

● G Gray

● B Blue




T



F	Fixed support VESA 75/100
R	Height adjustable support
N	Notebook enclosure
C	PC case with screw system
L	HPL work shelf
H	HPL work shelf with handles
M	Metal plate
20	200N max load 17 kg
30	300N max load 26 kg
40	400N max load 35 kg
T	TWIN WHEELS Ø 100 mm 4 with Brakes, 2 ESD
S	SINGLE WHEELS Ø 100 mm 4 with Brakes, 2 ESD

Platform	Color	Wheels	Column	Console	Monitor Support
X1T	G B	S T	20 30 40	N C L H M	/ F R
☐	☐ ☐	☐ ☐	☐ ☐ ☐	☐ ☐ ☐ ☐ ☐	☐ ☐ ☐



X1T Platform

General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	820 - 1170 mm
Center Vesa height range	1146 - 1496 mm
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.


Ref.	Color Handle Bumper	Wheel	Gas spring strenght	Upper column max load
X1T-BS-20L-F	Blue	Single 100mm	200N	17 kg
X1T-BT-20L-F	Blue	TWIN 100mm(2ESD)	200N	17 kg
X1T-BS-30L-F	Blue	Single 100mm	300N	26 kg
X1T-BT-30L-F	Blue	TWIN 100mm(2ESD)	300N	26 kg
X1T-BS-40L-F	Blue	Single 100mm	400N	35 kg
X1T-BT-40L-F	Blue	TWIN 100mm(2ESD)	400N	35 kg



 *AVAILABLE ALSO IN GREY



X1T Platform




General information	
Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with front handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	820 - 1170 mm
Center Vesa height range	1146 - 1496 mm
Vesa max load	12kg
Vesa swivel	(+/30)
Vesa tilting	-0°+40°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color Handle Bumper	Wheel	Gas spring strenght	Upper column max load
X1T-BS-20H-F	Blue	Single 100mm	200N	17 kg
X1T-BT-20H-F	Blue	TWIN 100mm(2ESD)	200N	17 kg
X1T-BS-30H-F	Blue	Single 100mm	300N	26 kg
X1T-BT-30H-F	Blue	TWIN 100mm(2ESD)	300N	26 kg
X1T-BS-40H-F	Blue	Single 100mm	400N	35 kg
X1T-BT-40H-F	Blue	TWIN 100mm(2ESD)	400N	35 kg

  *AVAILABLE ALSO IN GREY



X1T Platform


General information	
Base foot print	525x525mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	820 - 1170 mm
Center Vesa height range	1090 - 1570 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20L-R	Blue	Single 100mm	200N	17
X1T-BT-20L-R	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30L-R	Blue	Single 100mm	300N	26
X1T-BT-30L-R	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40L-R	Blue	Single 100mm	400N	35
X1T-BT-40L-R	Blue	TWIN 100mm(2ESD)	400N	35



*AVAILABLE ALSO IN GREY



X1T Platform

General information


Base foot print	520x520mm
Max Load	50kg
Top Shelf Size	510(W)X410(D)mm
Top Material	HPL with front handle
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	820 - 1170 mm
Center Vesa height range	1090 - 1570 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20H-R	Blue	Single 100mm	200N	17
X1T-BT-20H-R	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30H-R	Blue	Single 100mm	300N	26
X1T-BT-30H-R	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40H-R	Blue	Single 100mm	400N	35
X1T-BT-40H-R	Blue	TWIN 100mm(2ESD)	400N	35



*AVAILABLE ALSO IN GREY




X1T Platform

General information	
Base foot print	520x520mm
Case (S09) inner size	(W)440 (D)280 (H) 74 mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Top Material	HPL/PC case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	855 - 1205 mm
Center Vesa height range	1210 - 1560 mm
Vesa max load	12kg
Vesa swivel	(+/-30)
Vesa tilting	-0°+40°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20C-F	Blue	Single 100mm	200N	17
X1T-BT-20C-F	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30C-F	Blue	Single 100mm	300N	26
X1T-BT-30C-F	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40C-F	Blue	Single 100mm	400N	35
X1T-BT-40C-F	Blue	TWIN 100mm(2ESD)	400N	35

● ● *AVAILABLE ALSO IN GREY




X1T Platform

General information	
Base foot print	520x520mm
Case (S09) inner size	(W)440 (D)280 (H) 74 mm
Max Load	50kg
Top Shelf Size	(W)513 (D)406mm
Top Material	HPL/PC case
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	855 - 1205 mm
Center Vesa height range	1125 - 1605 mm
Vesa max load	12kg
Vesa swivel	(+/-90)
Vesa tilting	-15°+25°
Height adj	35 cm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20C-R	Blue	Single 100mm	200N	17
X1T-BT-20C-R	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30C-R	Blue	Single 100mm	300N	26
X1T-BT-30C-R	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40C-R	Blue	Single 100mm	400N	35
X1T-BT-40C-R	Blue	TWIN 100mm(2ESD)	400N	35

  *AVAILABLE ALSO IN GREY



X1T Platform

General information	
Base foot print	520x520mm
Max notebook dimension	17"
Max Load	50kg
Console Height range	856 - 1206 mm
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100mm wheel, the trolley is 20 mm shorter.


Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20N	Blue	Single 100mm	200N	17
X1T-BT-20N	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30N	Blue	Single 100mm	300N	26
X1T-BT-30N	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40N	Blue	Single 100mm	400N	35
X1T-BT-40N	Blue	TWIN 100mm(2ESD)	400N	35



 *AVAILABLE ALSO IN GREY



X1T Platform



General information	
Base foot print	525x525mm
Max Load	50kg
Top Shelf Size	440(W)X400(D)mm
Top Material	Steel
Compliant with	MDR 2017/745 EN 60601-1
Max Console Load	20kg
Console height range	812 - 1162 mm
Packaging	800x600x400 mm
Total Weight without Accessories	30 kg

When using the TWIN 100mm wheel, the trolley is 20 mm shorter.

Ref.	Accent Color	Wheel	Gas spring strenght	Upper part max load
X1T-BS-20M	Blue	Single 100mm	200N	17
X1T-BT-20M	Blue	TWIN 100mm(2ESD)	200N	17
X1T-BS-30M	Blue	Single 100mm	300N	26
X1T-BT-30M	Blue	TWIN 100mm(2ESD)	300N	26
X1T-BS-40M	Blue	Single 100mm	400N	35
X1T-BT-40M	Blue	TWIN 100mm(2ESD)	400N	35



*AVAILABLE ALSO IN GREY

all the accessories you need

Ask for custom accessories, we will
design them for your application!

IV POLE

ART. CODE: **015MLTACC**

(H) 350mm

Max load 4Kg
2Kg per hook



CAMERA HOLDER

ART. CODE: **010MLTACC**

**Support to be mounted on VESA with
Height adjustment from 19" to 27"**

Camera Plate 150x150 mm



Max load 1,5Kg



PIVOT HINGE FOR HEIGHT ADJ.VESA

ART. CODE: **019MLTACC**

**Support to be mounted on VESA
+-90°**




SIDE RAIL 10x25mm DIN rail

ART. CODE: **028MLTACC**

Length 30 cm

Max load 4Kg




PROBE HOLDER A

ART. CODE: **023MLTACC**

Inner ø min 16 mm
Inner ø max 55 mm
Inner Height 42mm

Max load 1Kg



PROBE HOLDER B

ART. CODE: **024MLTACC**

Inner ø min 20 mm
Inner ø max 55 mm
Inner Height 82mm

Max load 3Kg



SANITIZER GEL HOLDER

ART. CODE: **013MLTACC**

Inner dimensions:
80x80mm
Height 112mm


Max load 1Kg



BARCODE A

ART. CODE: **005MLTACC**


Max load 2Kg



BARCODE B

ART. CODE: **006MLTACC**

Max load 2Kg



PROBE HOLDER C

ART. CODE: **033MLTACC**

PMMA Transparent



Max load 3Kg

PROBE HOLDER D

ART. CODE: **034MLTACC**

PMMA Transparent



Max load 3Kg

PROBE HOLDER E

ART. CODE: **035MLTACC**

PMMA Transparent




Max load 3Kg

QWERTY SUPPORT

ART. CODE: **025MLTACC**

(W) 490mm
(D) 155mm

Max load 1Kg



L-R MOUSEPAD

ART. CODE: **017MLTACC**

REQUIRES QWERTY SUPPORT

(W) 140 mm
(D) 170 mm

Max load 1Kg



PC BASKET

ART. CODE: **007MLTACC**


Polycarbonate, finish: frosted



(W) 333mm
(D) 243mm
(H) 101mm
Max load 3Kg

WIRED BIN

ART. CODE: **031MLTACC**



(W) 315mm
(D) 100mm
(H) 135mm

Max load 3Kg

PC DRAWER

ART. CODE: **014MLTACC**

Polycarbonate, finish: frosted



(W) 333 mm
(D) 243mm
(H) 101mm
Max load 3Kg

SHELF

ART. CODE: **032MLTACC**




Internal dimensions
(W) 390mm
(D) 310mm

Max load 3Kg

COLUMN RAIL

ART. CODE: **011MLTACC**



(W) 300mm,
Max load 4Kg

PC SUPPORT

ART. CODE: **030MLTACC**



Internal dimensions
(W) 440mm
(D) 185mm

Max load 10Kg

PRINTER TRAY

ART. CODE: **021MLTACC**



Inner dim.
(W) 450 mm,
(D) 380mm,

Max load 10Kg

ISO DRAWER


ART. CODE: **004MLTACC**

with KEY LOCKING SYSTEM



Outer dim.
(W) 465 mm
(D) 360 mm
(H) 268mm

Max load 6Kg



removable tray ISO
40x30x10



MEDICAL SOCKET SUPPORT

ART. CODE: **018MLTACC**

**Ehman multi sockets not included
To be ordered separately
(medical multi sockets)**

(W) 96 mm
(H) 453 mm

Max load 3Kg



AL DRAWER K

ART. CODE: **003MLTACC**

with KEY LOCKING SYSTEM



Inner dim.
(W) 347mm
(D) 247mm
(H) 110 mm

Outer dim.
(W) 465 mm
(D) 360 mm
(H) 134mm

Max load 3Kg



with NUMERIC LOCKING SYSTEM




compatible with SALTO lockers

SOCKET/PSU SUPPORT

ART. CODE: **029MLTACC**

(W) 96 mm
(H) 238 mm

Max load 3Kg



DC POWER BATTERY SUPPORT


ART. CODE: **012MLTACC**

ready for QPD DC Box
DCPower4C battery

Housing

(W) 260 mm
(D) 105 mm
(H) 351mm


Max load 9Kg




Some examples of custom Configurations


Custom accessories for medical trolleys are essential for optimizing efficiency and organization in healthcare settings. These accessories, tailored to specific needs, include drawer dividers, lockable compartments, adjustable shelving, power outlets, waste bins, accessory hooks, documentation holders, and customized labeling. By integrating these features, medical facilities can streamline workflow, enhance patient care delivery, and ensure a clean and organized environment.

Customized console
for ultrasound system






Customized support
for endo devices
with inclination of 8°



Custom solutions




No console version
N. 2, 24" panel PC



Double monitor
N. 2, 24" panel PC





Custom solutions



MLTCAT X1T_03 (EN)



TouchPoint Group

Multix Solutions s.r.l.
phone +39 031 4100263

info@multix-trolley.com
www.multix-trolley.com

Via Padre Rovagnati, 2
22066 Mariano Comense
COMO - ITALY

Made in Italy



CATALOGUE

X1 T