



X2 Platform

ULTIMATE SIMPLICITY

MULTIX

X2 PLATFORM

Ultimate simplicity

Its new base with a central column ensures user comfort and swift manoeuvrability in fast-paced medical settings.

Multiple applications:


ideal for **Monitoring** but also for **ECG**, **Ultrasound**, **Ventilators** and **Point of Care**.

Agile: effortless mobility for enhanced workflow efficiency.

Customizable: Tailor dimensions and add specialized features to meet your specific needs.

The **X2 Platform** sets a new standard for efficiency and mobility in healthcare, transforming the way you deliver care.

Multix, Your Partner
in **Healthcare Excellence**





Innovative medical monitoring device mounting system

Quick and Easy Installation:

Features a spring-loaded snap mechanism for fast and secure attachment.

Effortless Removal:

Devices can be detached simply by pulling a lever.

Safe and Reliable:


Ensures a secure hold, ideal for medical environments.

User-Friendly Design:


Practical, intuitive, and convenient for monitoring vital signals

Ergonomics:

Tilting $-10^{\circ} + 20^{\circ}$
Double handle



Adaptor plates for Monitoring



Let us inspire you

X2 PLATFORM


Modularity

The modular design approach enables customization and scalability, allowing the system to be configured for a variety of applications, including ECG, ventilators, and more.

Each module can be easily adjusted or replaced, offering flexibility to adapt the cart to evolving medical needs.

This ensures versatility, efficiency, and seamless integration across different healthcare environments.


Multix, Your Partner
in **Healthcare Excellence**



Accent Color

● **G** Gray

● **B** Blue



I	IV Pole
T	Tablet support
F	Vesa pole
L	Column 1100mm
M	ECG/device console
H	HPL console
V	Monitoring console
P	Handle
70	Column 694mm
84	Column 838mm
110	Column 1100mm
X2	Base Aluminium 520x520mm

[illegible]

X2

Specifications

Following are the main specifications of the system, highlighting its key features and performance standards. These specifications ensure reliability, adaptability, and ease of use in a wide range of medical applications.

Compliant with:
MDR 2017/745
EN 60601-1

Multix, Your Partner
in Healthcare Excellence


X2 Easy Fix

General information	
Base foot print	520x520mm
Max Load	44kg
Max Console Load	14kg
Console Height	1120 mm
Vesa tilting	-10°+20°
Total Weight without Accessories	9.8kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color	Wheel ø100mm	Column Lenght	Console Height	Packaging Dim. mm
X2N-BS-84V	Blue	Single	838mm	1040 mm	880x600x370
X2N-BT-84V	Blue	Twin	838mm	1040 mm	880x600x370

  AVAILABLE ALSO IN GREY



X2 standard metal console

General information	
Base foot print	520x520mm
Max Load	44kg
Max Console Load	25kg
Console Height	1031 mm
Console Size	380x380mm
Total Weight without Accessories	12.6 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color	Wheel ø100mm	Column Lenght	Console Height	Packaging Dim. mm
X2N-BS-84M	Blue	Single	838mm	1040 mm	880x600x370
X2N-BT-84M	Blue	Twin	838mm	1040 mm	880x600x370

  AVAILABLE ALSO IN GREY



X2 standard HPL console



General information	
Base foot print	520x520mm
Max Load	44kg
Max Console Load	25kg
Console Height	1040 mm
Console Size	500 x 480 mm
Total Weight without Accessories	10.2 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color	Wheel ø100mm	Column Lenght	Console Height	Packaging Dim. mm
X2N-BS-84H	Blue	Single	838mm	1040 mm	880x600x370
X2N-BT-84H	Blue	Twin	838mm	1040 mm	880x600x370

  AVAILABLE ALSO IN GREY




X2-Monitor

General information	
Base foot print	520x520mm
Max Load	44kg
Compliant with	MDR 2017/745 EN 60601-1
Vesa max load	14kg
Vesa tilting	-0°+35°
Vesa swivel	(+/-28°)
Total Weight without Accessories Includes 2 Counterweight.	13.7 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color	Wheel ø100mm	Column Lenght	C. VESA Height	Packaging Dim. mm
X2-BS-84-P-F	Blue	Single	838mm	1510 mm	880x600x370
X2-BT-84-P-F	Blue	Twin	838mm	1510 mm	880x600x370
X2-BS-110-P-F	Blue	Single	1100mm	1770 mm	1200x800x280
X2-BT-110-P-F	Blue	Twin	1100mm	1770 mm	1200x800x280

  AVAILABLE ALSO IN GREY



X2-Tablet

General information	
Base foot print	520x520mm
Max Load	44kg
Compliant with	MDR 2017/745 EN 60601-1
Vesa Max Load	6 kg
Vesa Tilting	-0°+35°
Vesa swivel	(+/-28°)
Tablet size	9" up to 11"
Total Weight without Accessories	10.3 kg

When using the TWIN 100 mm wheel, the trolley is 20 mm lower.

Ref.	Color	Wheel ø100mm	Column Lenght	C. VESA Height	Packaging Dim. mm
X2N-BS-84-P-T	Blue	Single	838mm	1370 mm	880x600x370
X2N-BT-84-P-T	Blue	Twin	838mm	1370 mm	880x600x370
X2N-BS-110-P-T	Blue	Single	1100mm	1631mm	1200x800x280
X2N-BT-110-P-T	Blue	Twin	1100mm	1631mm	1200x800x280

  AVAILABLE ALSO IN GREY

Ref.	Color	Wheel ø100mm	Column Lenght	Packaging Dim. mm
X2-BS-70-P-I	Blue	Single	694mm	800x570x370
X2-BT-70-P-I	Blue	Twin	694mm	800x570x370
X2-BS-84-P-I	Blue	Single	838mm	880x600x370
X2-BT-84-P-I	Blue	Twin	838mm	880x600x370

  AVAILABLE ALSO IN GREY

X2-Infusion



General information	
Base foot print	520x520mm
Max Load	44kg
Compliant with	MDR 2017/745 EN 60601-1
Max IV Pole Load	2 kg per hook
IV Pole Height range	1380-1900 mm
Total Weight without Accessories Includes 2 Counterweight.	13.5 kg


When using the TWIN 100 mm wheel, the trolley is 20 mm lower.



Some examples of standard and custom Configurations

Custom accessories for medical trolleys are essential for optimizing efficiency and organization in healthcare settings. These accessories, tailored to specific needs, include adjustable shelving, power outlets, waste bins, accessory hooks and customized labeling. By integrating these features, medical facilities can streamline workflow, enhance patient care delivery, and ensure a clean and organized environment.

Multix, Your Partner
in **Healthcare Excellence**





Custom Solution

based on X2 platform
Custom console and holder
for specific device

Wired Bin

(031MLTACC)

Double 10L Cylinder Support

(076MLTACC)

Custom Solution

based on X2 platform
Custom console, handle for
NKV 330 Ventilator

Wired Basket With Cable holder

(054MLTACC)

Double 5L Cylinder Support

(077MLTACC)

Humidifier support

(079MLTACC)



Let us inspire you

MULTIX



Scan the brochure



TOUCH POINT GROUP

Multix Solutions s.r.l.
phone +39 031 4100263

info@multix-trolley.com
www.multix-trolley.com

Via Padre Rovagnati, 2
22066 Mariano Comense
COMO - ITALY

Made in Italy



X2 Platform
CATALOGUE