





Overall dimensions and weight:	550 mm (W) × 781 mm (D) × 1265 mm (H min) - 45 kg
Work Height	Adjustable with electrical column system 800 mm (H min) - 1220 mm (H max) - 420mm stroke
Console	Polycarbonate plastic RAL7047 - 10 kg max total load (top+internal)
Top console	Color White 410 with integrated front handles
Top console dimension	513 mm (W) × 495 mm (D) x 12 mm (H) -
Internal technical compartment dimension	450 mm (W) × 180 mm (D) x 74 mm (H)
With integrated internal medical grade power strip	5x DIN C13 power outlet-1x DIN C14 power inlet-6x Potential equalization node Max rating 250 Vac - 10A (250Vac - 15A UL&CSA) Main power cable NOT included. Available on request in the format of the country of use Illuminated main switch on the back of the console
With ON/OFF button and battery charge status indicator for the height adjustment system	With 4 green LED of charge status (25%-50%-75%-100%)
Handle	Die-cast aluminum overmolded with silicone Blue 285
Pole	Painted steel - Dark grey RAL7016 - 12 kg max load
System status indicator LED	Green Light
Console height adjustment control	Circular up and down button
VESA 75-100 attachment	With tilting (0°/+30°) and swivelling (-90°/+90°)
Column	Painted aluminium - Grey RAL7047
Base	Steel with polypropylene plastic covers with Blue 285 silicone bumper
Base footprint	548 mm (W) × 578 mm (D)
4 twin wheels with brake	ø 125 mm (2x ESD)





Optional accessories

003MLTACC - Aluminium drawer (compatible SALTO)
004MLTACC - Aluminium ISO drawer
005MLTACC - Barcode
006MLTACC - Barcode 2
007MLTACC - Basket
008MLTACC - Cable channel 470
010MLTACC - Camera holder
011MLTACC - Column rail
013MLTACC - Disinfectant gel holder
014MLTACC - Drawer
015MLTACC - IV pole
017MLTACC - L-R mousepad
022MLTACC - Printer tray
023MLTACC - Probe holder
024MLTACC - Probe holder 2
027MLTACC - Sharp container holder
028MLTACC - Side rail
029MLTACC - Socket support
030MLTACC - Tower PC support
031MLTACC - Wired bin
032MLTACC - Universal tray

	033MLTACC - Multiple probe holder d.30-d.22
	034MLTACC - Multiple probe holder d.22-d.15
	035MLTACC - Multiple probe holder d.30-d.15
	067MLTACC - Shelf 8°
	073MLTACC - Multisocket 4 DIN
	075MLTACC - Plastic drawer
	080MLTACC - X3E QWERTY support
Total weight included safe working load	67 Kg
Packaging dimensions and weight	757 mm (W) × 661 mm (D) × 1100 mm (H) - 65 kg
CE Declaration of Conformity	EN 60601-1-2 2017/745 MDR





Overall dimensions and weight:	550 mm (W) × 761 mm (D) × 810 mm (H min) - 42 kg
Work Height	Adjustable with electrical column system 800 mm (H min) - 1220 mm (H max) - 420mm stroke
Console	Polycarbonate plastic RAL7047 - 10 kg max total load (top+internal)
Top console	Color White 410
Top console dimension	513 mm (W) × 405 mm (D) x 12 mm (H)
Internal notebook compartment dimension	For notebook from 13" to 17"
With ON/OFF button and battery charge status indicator for the height adjustment system	With 4 green LED of charge status (25%-50%-75%-100%)
Handle	Die-cast aluminum overmolded with silicone Blue 285
Top column control panel	Nylon plastic black
System status indicator LED	Green Light
Console height adjustment control	Circular up and down button
Column	Painted aluminium - Grey RAL7047
Base	Steel with polypropylene plastic covers with Blue 285 silicone bumper
Base footprint	548 mm (W) × 578 mm (D)
4 twin wheels with brake	ø 125 mm (2x ESD)
Keyboard support	Painted aluminium - Color dark grey 432 Painted steel - Color dark grey 432
Keyboard max dimensions	450 mm (W) x 133 mm (D) x 35-50mm (H) - 1 kg max load



Optional accessories

003MLTACC - Aluminium drawer (compatible SALTO)
004MLTACC - Aluminium ISO drawer
005MLTACC - Barcode
006MLTACC - Barcode 2
007MLTACC - Basket
008MLTACC - Cable channel 470
010MLTACC - Camera holder
011MLTACC - Column rail
013MLTACC - Disinfectant gel holder
014MLTACC - Drawer
015MLTACC - IV pole
017MLTACC - L-R mousepad
018MLTACC - Medical socket support
022MLTACC - Printer tray
023MLTACC - Probe holder
024MLTACC - Probe holder 2
026MLTACC - Reomed 600 support
027MLTACC - Sharp container holder
028MLTACC - Side rail
029MLTACC - Socket support
030MLTACC - Tower PC support
031MLTACC - Wired bin
032MLTACC - Universal tray
033MLTACC - Multiple probe holder d.30-d.22
034MLTACC - Multiple probe holder d.22-d.15

035MLTACC - Multiple probe holder d.30-d.15

067MLTACC - Shelf 8°

073MLTACC - Multisocket 4 DIN

075MLTACC - Plastic drawer

Total weight included safe working load

52 Kg

Packaging dimensions and weight

757 mm (W) × 661 mm (D) × 1100 mm (H) - 62 kg


CE Declaration of Conformity

EN 60601-1-2

2017/745 MDR



Overall dimensions and weight:	550 mm (W) × 781 mm (D) × 1265 mm (H min) - 52 kg (with 720Wh battery)
Work Height	Adjustable with electrical column system 800 mm (H min) - 1220 mm (H max) - 420mm stroke
Console	Polycarbonate plastic RAL7047 - 10 kg max total load (top+internal)
Top console	Color White 410 with integrated front handles
Top console dimension	513 mm (W) × 495 mm (D) x 12 mm (H)
Internal technical compartment dimension	450 mm (W) × 180 mm (D) x 74 mm (H)
With integrated internal medical grade power strip	5x DIN C13 power outlet-1x DIN C14 power inlet-6x Potential equalization node Max rating 250 Vac - 10A (250Vac - 15A UL&CSA) directly connected to the DC Power battery Illuminated main switch on the back of the console
With ON/OFF button and battery charge status indicator for the height adjustment system	With 4 green LED of charge status (25%-50%-75%-100%)
With integrated LED lamp to illuminate the keyboard or the floor	220mm 4000K white LED lamp with proximity ON/OFF sensor switch
Handle	Die-cast aluminum overmolded with silicone Blue 285
Pole	Painted steel - Dark grey RAL7016 - 12 kg max load
Status indicator LED	Green Light
Console height adjustment control	Circular up and down button
DC Power battery status display integrated support	2.8" Graphic display with full state of charge and remaining runtime
VESA 75-100 attachment	With tilting (0°/+30°) and swivelling (-90°/+90°)



Column	Painted aluminium - Grey RAL7047
Base	Steel with polypropylene plastic covers with Blue 285 silicone bumper
Base footprint	548 mm (W) × 578 mm (D)
4 twin wheels with brake	Ø 125 mm (2x ESD)
DC Power battery dedicated support	DC Power QPD BOX support
DC Power battery available versions	240Wh battery-2 DC/DC converter with internal inverter for 110-230 Vac output -4.2Kg 480Wh battery-2 DC/DC converter with internal inverter for 110-230 Vac output -5.7Kg 720Wh battery-2 DC/DC converter with internal inverter for 110-230 Vac output -7.1Kg
Optional accessories	003MLTACC - Aluminium drawer (compatible SALTO) 004MLTACC - Aluminium ISO drawer 005MLTACC - Barcode 006MLTACC - Barcode 2 007MLTACC - Basket 008MLTACC - Cable channel 470 010MLTACC - Camera holder 011MLTACC - Column rail 013MLTACC - Disinfectant gel holder 014MLTACC - Drawer 015MLTACC - IV pole 017MLTACC - L-R mousepad 018MLTACC - Medical socket support

	022MLTACC - Printer tray
	023MLTACC - Probe holder
	024MLTACC - Probe holder 2
	025MLTACC - Keyboard support
	026MLTACC - Reomed 600 support
	027MLTACC - Sharp container holder
	028MLTACC - Side rail
	029MLTACC - Socket support
	030MLTACC - Tower PC support
	031MLTACC - Wired bin
	032MLTACC - Universal tray
	033MLTACC - Multiple probe holder d.30-d.22
	034MLTACC - Multiple probe holder d.22-d.15
	035MLTACC - Multiple probe holder d.30-d.15
	067MLTACC - Shelf 8°
	073MLTACC - Multisocket 4 DIN
	075MLTACC - Plastic drawer
	080MLTACC - X3E QWERTY support
Total weight included safe working load	74 Kg
Packaging dimensions and weight	757 mm (W) × 661 mm (D) × 1100 mm (H) - 65 kg
CE Declaration of Conformity	EN 60601-1-2 2017/745 MDR

