



Overall dimensions and weight:	550 mm (W) × 642 mm (D) × 1307 mm (H min) - 34 kg
Work Height	Adjustable with mechanical gas spring system 891 mm (H min) - 1136 mm (H max)
Top Console in HPL Laminate	Color White ABET PRINT 410
Handle in Silicone-overmolded aluminum	Color Blue PANTONE 285
Top Console dimension	513 mm (W) × 406 mm (D) x 12 mm (H) - 10 kg max load
Pole with VESA 75-100 attachment	With tilting (0°/+38°) and swivelling (-30°/+30°) - 8 kg max load
Aluminium column	with lateral one-way cable channels
Steel base with PP cover	with Blue PANTONE 285 silicone bumpers
Base footprint	548 mm (W) × 578 mm (D)
4 twin wheels with brake	ø 125 mm (2ESD - 2STD)
Optional accessories	
	001MLTACC - Advantech battery support X3
	002MLTACC - Aluminium drawer
	003MLTACC - Aluminium drawer (compatible SALTO)
	004MLTACC - Aluminium ISO drawer
	005MLTACC - Barcode
	006MLTACC - Barcode 2
	007MLTACC - Basket

	008MLTACC - Cable channel 470
	010MLTACC - Camera holder
	011MLTACC - Column rail
	012MLTACC - DC Power battery support
	013MLTACC - Disinfectant gel holder
	014MLTACC - Drawer
	015MLTACC - IV pole
	017MLTACC - L-R mousepad
	018MLTACC - Medical socket support
	019MLTACC - Pivot hinge for height adjustable vesa pole
	022MLTACC - Printer tray X3-X3E
	023MLTACC - Probe holder
	024MLTACC - Probe holder 2
	025MLTACC - Keyboard support
	026MLTACC - Reomed 600 support X3
	027MLTACC - Sharp container holder
	028MLTACC - Side rail
	029MLTACC - Socket support
	030MLTACC - Tower PC support
	031MLTACC - Wired bin
	032MLTACC - Universal tray
	033MLTACC - Multiple probe holder d.30-d.22
	034MLTACC - Multiple probe holder d.22-d.15
	035MLTACC - Multiple probe holder d.30-d.15
Total weight	80 Kg
Packaging dimensions and weight	757 mm (W) × 661 mm (D) × 985 mm (H) - 70 kg
CE Declaration of Conformity	EN 60601-1 2017/745 MDR