



### X1-BT-77C-F

General information	
<b>Base foot print</b>	520x520mm
<b>Max Load</b>	50kg
<b>Top Shelf Size</b>	(W)513 (D)406mm
<b>Case (S09) inner size</b>	(W)440,(D)280,(H)74 mm
<b>Top Material</b>	HPL with plastic case
<b>Compliant with</b>	MDR 2017/745 EN 60601-1
<b>Max Console Load</b>	20kg
<b>Vesa max load</b>	12kg
<b>Vesa swivel</b>	(+/-30)
<b>Vesa tilting</b>	-0°+40°

\*AVAILABLE ALSO IN GREY